



**SYMMETRIC
FUMONISIN GREEN**

LATERAL FLOW TEST KIT

for the quantitative determination of Fumonisin in grains, cereals and animal feed

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric Fumonisin Green, S7024/S7048, is a Lateral Flow Test kit for the quantitative determination of Fumonisin in grains, cereals and animal feed.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Corn, Corn flour, Wheat, Wheat flour, Barley, Malt, Oats, Soy beans, Soybean meal, DDGS, Sunflower Meal, Brown rice, White rice, Rice Flour, Buckwheat, Millet, Dried Brassica Integrifolia, Dried Gai Choy, Dried Palm, Pasta, Pop corn, Beer residue

- Sample preparation: extraction
- Test time (reaction time after samples and reagents preparation): 3min
- Range: 0 - 4ppm
- Shelf life: 12 months
- Storage: 2-8°C

Specifications

- The LOD of the method is 0.1ppm Fumonisin.
- The LOQ of the method is 0.15ppm Fumonisin.
- Cross-reactivity: The cross-reaction of the anti-Fumonisin antibody with FB1, FB2 and FB3 is 100, 65 and 48% respectively.

1. Description

Symmetric Fumonisin Green is an innovative Lateral Flow test, utilizing state-of-the-art features for the quantitative detection of Fumonisin in grains, cereals and animal feed. This Lateral Flow test utilizes an ecological solution [1-5] for the extraction step, instead of the usual organic solvents.

2. General Information

Fumonisin is a member of the trichothecene mycotoxins produced by fungi of *Fusarium moniliforme* (*F. verticillioides*), *F. proliferatum*, and several other *Fusarium* species. Grains including corn, wheat and other cereals are frequently infected by these fungi in the field or during storage. More than ten types of fumonisins have been isolated and characterized. Of these, Fumonisin B1 (FB1), B2 (FB2), and B3 (FB3) are the major fumonisins produced. Fumonisin is hepatotoxic and nephrotoxic in all animal species tested. The earliest histological change to appear in either the liver or kidney of fumonisin-treated animals is increased apoptosis followed by regenerative cell proliferation, while the acute toxicity of fumonisin is low, it is the known cause of two diseases which occur in domestic animals with rapid onset: equine leukoencephalomalacia and porcine pulmonary oedema syndrome. Both of these diseases involve disturbed sphingolipid metabolism and cardiovascular dysfunction. Most controlling government agencies worldwide have regulations regarding the amount of FB1, and FB2 allowable in human and animal foodstuffs. Accurate and rapid determination of Fumonisin presence in commodities is of paramount importance.

3. Principle of the Method

The quantitative lateral flow test is based on the immunochromatography assay principles. The wells of the microtiter strips contain Fumonisin specific antibodies conjugated to colloidal gold. Diluted extract is added into the well. A dipstick with two capture lines, test and control, is dipped into the well. The suspended mixture starts flowing vertically on the dipstick and passes through the two lines. While running, Fumonisin (if it is present) binds to the antibodies. A valid test should always have the upper control line red. If the sample is free of Fumonisin, a color development occurs at the test line, indicating the absence of Fumonisin in the sample. On the contrary, the presence of Fumonisin in the sample will cause a reduced colored signal at the test line. The test line color intensity is indirectly proportionate to the concentration of Fumonisin present in the samples. By utilizing S-Flow software and the symmetric quantification technology [6, 7], Fumonisin is accurately quantified.

4. Reagents Provided

Symmetric Fumonisin Green kit contains sufficient reagents and materials for 24/48 reactions.

Reagents (Store at 2-8°C)	Quantity for 24 reactions	Quantity for 48 reactions
Pots each with 1 strip of 8 reagent microwells and 8 dipsticks	3	6
Sample Diluent Tubes	24	48
Extraction Solution 10X (50ml)	1	2
High Range Solution (10ml)	1	1

5. Materials required but not provided

- A grinder sufficient to render sample to particle size of fine instant coffee
- Balance with 0 - 50g measuring capability and Graduated cylinder - 50ml
- Deionized water
- Mini centrifuge (spin) and plastic tubes 1,5 or 2ml
- Tube roller or Vortex mixer
- 100 or 200µl adjustable micropipettes (single or multi channel) with disposable tips
- S-Flow** software along with matching scanner device

6. Storage Instructions

Store kit components between 2 - 8°C. Do not freeze any components provided. Reseal the unused strips in the storing tube together with the desiccant bag provided. The expiry date of the kit and reagents is stated on their labels and no quality guarantee is accepted after the expiration date. The expiry of the kit components can only be guaranteed if the components are stored properly and the reagent is not contaminated due to prior handling. Do not interchange individual components between kits of different lot numbers.

7. Safety and Precautions for use

Let the reagents warm to room temperature (21 - 25°C) before the analysis (at least half an hour) and cover them when not in use. Use a clean disposable plastic pipette tip for each reagent, to avoid cross contamination.

8. Preparation of Extraction Solution

In case of the occurrence of crystals in the **Extraction Solution 10X**, the warming by gentle dismantling (using hands) of the crystals is needed. Pour entire content of the solution concentrate (50ml) into a clean 500ml graduated cylinder, rinse the vial with distilled or deionized water and pour the content again into the cylinder and fill to a final volume of 500ml with distilled or deionized water (50ml Extraction Solution 10X and 450ml deionized water). Mix gently to avoid foaming, transferring the final solution from cylinder to a clean bottle and back two times. The clean bottle with **1X Extraction Solution** working solution can be left out of the refrigerator during the method procedure and subsequent be stored 2 - 8°C for one year.

9. Sample preparation

- The sample must be collected according to established sampling techniques. Grind a representative sample to the particle size of fine instant coffee (50% passes through a 20 mesh screen).
- Weigh out a 5g ground portion of the sample and add 15ml of the Extraction Solution (see 8). Mix using a tube roller for 5 minutes (or vortex for 2min). **The ratio of sample to Extraction Solution is 1:3 (w/v).**
- Allow the particulate matter to settle. Centrifuge 1ml of the extract for 2min using a mini centrifuge (spin). (The extracted sample should have pH value of 6.2 - 7.0. If the pH is less than 6.2, the pH should be neutralized using NaOH.)
- Add **100µl** of extract (supernatant) into the Sample diluent tube provided and mix well. Run the diluted extract within **30 minutes**.

Note: In case the result is greater than 4000 ppb, the sample should be further diluted with **High Range Solution** and re-tested. To achieve a dilution factor of 5 or 10, add 100µl of the already diluted sample into 400µl or 900µl of **High Range Solution (respectively)**. Run the second dilution within **30 minutes**.

Set the suitable dilution factor type to multiply the results by **5** or **10**.

DILUTION FACTOR 5 (Quantify from 0,75 up to 20ppm)
100µl of the already diluted sample + 400µl of High Range Solution

DILUTION FACTOR 10 (Quantify from 1.5up to 40ppm)
100µl of the already diluted sample + 900µl of High Range Solution

10. Method Procedure

- Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the ambient temperature.
- Download and/or set the kit's **lot number**, as provided in the Quality Assurance Certificate and then set the suitable **Dilution Factor type**.
- Open one plastic pot and take out as many test strips and microwells as samples to be tested.
- The pot with dipsticks should **always be well closed** after reagents have been taken out.
- Dispense **100µl of diluted extract** into the microwell and pipette **up and down 4 times** to completely mix the lyophilized gold particles in the sample, while avoiding bubbles. The sample should turn into a **uniform pink color**. In case of more than 2 samples, an 8 channel multipipette should be used.
- Place the appropriate number of sticks into microwells **immediately** and set timer for 3 minutes.
- When the 3 minutes are over, take the dipsticks out of the microwells and remove the white cotton sample-pad of the stick **immediately**. Touch the stick with your hand from the colorful pad and remove the white pad with your hands. Do not use a paper towel or any other material.
- Place the stick inside the plastic holder in order to be scanned. In case of EPSON scanner, the **sticks must be facing down (inverted)** and the colored side must be facing the orange sticker. **NOTE:** The sticks should be scanned within 2 minutes after the sample-pads removal.
- The software will use a Lot specific curve to calculate the results (ppb). A simple visual interpretation of the stick is NOT possible.

11. Performance Evaluation

11.1 Reference Materials

Several reference materials are being used for the evaluation of each product of ProGnosis Biotech S.A. in the context of Quality Control performed by Quality Control Department. Please request a validation report, including the results, at info@prognosis-biotech.com.

11.2 Proficiency Tests

All products participate frequently in Proficiency Tests. For more information, visit the individual product page in our website: www.prognosis-biotech.com.

12. Method Summary

Total method time: 3 minutes

Extract the samples



Add 100µl of extract (supernatant) into the Sample diluent tube provided



Dispense 100µl of each sample into the microwells and mix 5 times the sample with the lyophilized gold particles



Place the appropriate number of sticks into microwells immediately.



(Wait 3 mins)

Take the stick out and remove the white sample-pad immediately



Place the stick in the appropriate device to be scanned



Quantify through s-flow software

VERSION N2

CAT.NUMBER: S7024/S7048

STORAGE: 2-8°C



LATERAL FLOW TEST KIT

for the quantitative determination of Fumonisin in grains, cereals and animal feed

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric Fumonisin Green, S7024/S7048, is a Lateral Flow Test kit for the quantitative determination of Fumonisin in grains, cereals and animal feed.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Corn, Corn flour, Wheat, Wheat flour, Barley, Malt, Oats, Soy beans, Soybean meal, DDGS, Sunflower Meal, Brown rice, White rice, Rice Flour, Buckwheat, Millet, Dried Brassica Integri-fovia, Dried Gai Choy, Dried Palm, Pasta, Pop corn, Beer residue

- Sample preparation: extraction
- Test time (reaction time after samples and reagents preparation): 3min
- Range: 0 - 4ppm
- Shelf life: 12 months
- Storage: 2-8°C

All immune assays supplied by ProGnosis Biotech S.A., are warranted to meet or exceed our published specification when used under normal conditions in your laboratory. If the product fails during the stated period, a replacement product will be issued.

ProGnosis Biotech S.A. makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. ProGnosis Biotech S.A. shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product. This method is considered to be a screening method, before a legal action, samples detected as positives must be confirmed with a confirmation method. This product is meant to be used only For Research or Manufacturing use and by qualified technicians.



www.prognosis-biotech.com
e: info@prognosis-biotech.com
t: +30 2410 623922
f: +30 700 700 6262
Farsalon 153 | 41335 Larissa, Greece

