


# VALIDATION REPORT

## SYMMETRIC FUMONISIN GREEN



Product code: S7024/S7048  
Document no: S70[V12]  
Date: 09 July 2024



## Symmetric FUMONISIN Lateral Flow kit

### INDEX

	<b>Page</b>
<b><u>1. Immunoassay Specifications</u></b>	
1.1 General Specifications	3
1.2 Specificity & Cross-reactivity	3
<b><u>2. Validation</u></b>	
2.1 Determination of the Limit of Detection LOD and the Limit of Quantification LOQ	3
2.2 Determination of Recovery (%)	3
<i>i. Determination of Recovery (%) at the LOQ level</i>	3
<i>ii. Determination of Recovery (%) for all matrices at two different levels</i>	4
2.3 Reproducibility	13
2.4 Performance Evaluation	13
<b><u>3. References</u></b>	13

## Symmetric FUMONISIN Lateral Flow kit

### 1. Immunoassay Specifications

#### 2.1 General Specifications

- The LOD of the method is 100 ppb (Table 1)
- The LOQ of the method is 150 ppb (Table 1)
- **Matrices:** Barley, Beer Residue, Brown Rice, Buckwheat, Corn, Corn Flour, DDGS, Dried Brassica Integrifolia, Dried Gai Choy, Dried Palm, Malt, Millet, Oats, Pasta, Pop Corn, Rice Flour, Soybeans, Soy bean meal, Sunflower meal, Wheat, Wheat Flour, White Rice.

#### 1.2 Specificity & Cross-reactivity

The cross-reaction of the anti-Fumonisin antibody with FB1, FB2 and FB3 is 100, 65 and 48% respectively.

### 2. Validation

#### 2.1 Determination of the Limit of Detection LOD and the Limit of Quantification LOQ

For the determination of LOD (2xSD) and LOQ (3xSD), two Fumonisin-free maize powder samples (<150 ppb) were used (Table 1).

**Table 1. Fumonisin-free maize powder samples for the determination of LOD and LOQ**

Sample (n=8)	Concentration (ppb)	
	MEAN	SD
Maize Powder A	7.16	49.70
Maize Powder B	5.75	49.81
	<b>SUM</b>	
MEAN	<b>6.46</b>	
SD	<b>49.76</b>	

The LOD and LOQ were defined as 2 x Standard Deviation and 3 x Standard Deviation of the Fumonisin-free maize powder blank samples, respectively. It was found that calculated LOD and LOQ is 100 ppb and 150 ppb, respectively.

#### 2.2 Determination of Recovery (%)

##### *Spike Protocol*

All the samples were spiked according to the United States Department of Agriculture (USDA). Agricultural Marketing Service. GIPSA's Federal Grain Inspection Service (FGIS) protocol. More particularly, all spike experiments were carried out by spiking the individual pre-weighted test portion with a concentration adjusted solution to maintain the spiking volume at 100  $\mu$ L. Liquid spike prepared in 100% methanol was added with a positive displacement syringe and dried for 30 minutes at 37°C prior to extraction, unless stated otherwise. A known amount of Aflatoxin Fumonisin is added to the solid sample to be tested using a standard solution. The dried spiked sample was then extracted and analyzed according to the manual S7024/S7048 V23 chapter 9.

##### *i. Determination of Recovery (%) at the LOQ level*

For the determination of Recovery at LOQ level (Table 2) Fumonisin-free maize powder B was used as blank and it has been spiked with Fumonisin Trilogy standard solution (TSL-204).

**Table 2. Recovery at LOQ level. Maize Powder B was spiked with Fumonisin.**

Sample (n=8)	Concentration (ppb)	Spike (ppb)	Recovery (%)	
MAIZE POWDER B	145.92	150.00	97.28	Average 150 ppb
			15.2	SD
			15.6	CV%

ii. Determination of Recovery (%) for all matrices at two different levels

For the determination of Recovery at two different levels (450 and 1500 ppb), different Fumonisin-free matrices were spiked with Fumonisin Trilogy standard solution (TSL-204) (Tables 3-24).

**Table 3. Recovery of Barley sample at two different levels.**

Spike 450 ppb		
Barley (n=8)	438.75	Average
	9.93	SD
	2.26	CV(%)
	97.5	Recovery (%)
Spike 1500 ppb		
Barley (n=8)	1457.31	Average
	45.12	SD
	3.09	CV(%)
	97.15	Recovery (%)
	97.32	Mean Recovery (%)

**Table 4. Recovery of Beer Residue sample at two different levels.**

Spike 450 ppb		
Beer Residue (n=8)	463.75	Average
	5.16	SD
	1.11	CV(%)
	103.05	Recovery (%)
Spike 1500 ppb		
Beer Residue (n=8)	1492.12	Average
	43.1	SD
	2.88	CV(%)
	99.47	Recovery (%)
	101.26	Mean Recovery (%)

**Table 5. Recovery of Brown Rice sample at two different levels.**

Spike 450 ppb		
Brown Rice (n=8)	463.25	Average
	8.88	SD
	1.9	CV(%)
	102.9	Recovery (%)
Spike 1500 ppb		
Brown Rice (n=8)	1461	Average
	67.22	SD
	4.6	CV(%)
	97.4	Recovery (%)
	100.15	Mean Recovery (%)

**Table 6. Recovery of Buckwheat sample at two different levels.**

Spike 450 ppb		
Buckwheat (n=8)	461.3	Average
	10.23	SD
	2.21	CV(%)
	102.51	Recovery (%)
Spike 1500 ppb		
Buckwheat (n=8)	1463.5	Average
	79.21	SD
	3.22	CV(%)
	97.56	Recovery (%)
	100.03	Mean Recovery (%)

**Table 7. Recovery of Corn sample at two different levels.**

Spike 450 ppb		
Corn (n=8)	461.25	Average
	25.52	SD
	5.53	CV(%)
	102.5	Recovery (%)
Spike 1500 ppb		
Corn (n=8)	1479.51	Average
	75.81	SD
	5.1	CV(%)
	98.63	Recovery (%)
	100.56	Mean Recovery (%)

**Table 8. Recovery of Corn flour sample at two different levels.**

Spike 450 ppb		
Corn flour (n=8)	481	Average
	9.08	SD
	1.88	CV(%)
	106.88	Recovery (%)
Spike 1500 ppb		
Corn flour (n=8)	1565.2	Average
	75.21	SD
	4.8	CV(%)
	104.3	Recovery (%)
	105.59	Mean Recovery (%)

**Table 9. Recovery of DDGS sample at two different levels.**

Spike 450 ppb		
Corn flour (n=8)	451.75	Average
	6.68	SD
	1.48	CV(%)
	100.4	Recovery (%)
Spike 1500 ppb		
Corn flour (n=8)	1492.12	Average
	57.11	SD
	3.82	CV(%)
	99.47	Recovery (%)
	99.93	Mean Recovery (%)

**Table 10. Recovery of Dried Brassica Intergrifolia sample at two different levels.**

Spike 450 ppb		
Dried Brassica Intergrifolia (n=8)	383.87	Average
	20.23	SD
	5.27	CV(%)
	85.3	Recovery (%)
Spike 1500 ppb		
Dried Brassica Intergrifolia (n=8)	1439.5	Average
	65.83	SD
	4.57	CV(%)
	95.96	Recovery (%)
	90.63	Mean Recovery (%)

**Table 11. Recovery of Dried Gai Choy sample at two different levels.**

Spike 450 ppb		
Dried Gai Choy (n=8)	386.12	Average
	20.24	SD
	5.24	CV(%)
	85.80	Recovery (%)
Spike 1500 ppb		
Dried Gai Choy (n=8)	1475.2	Average
	54.15	SD
	3.67	CV(%)
	97.14	Recovery (%)
	91.47	Mean Recovery (%)

**Table 12. Recovery of Dried Palm sample at two different levels.**

Spike 450 ppb		
Dried Palm (n=8)	388.12	Average
	19.2	SD
	4.95	CV(%)
	86.25	Recovery (%)
Spike 1500 ppb		
Dried Palm (n=8)	1448.33	Average
	57.1	SD
	3.94	CV(%)
	96.55	Recovery (%)
	91.4	Mean Recovery (%)

**Table 13. Recovery of Malt sample at two different levels.**

Spike 450 ppb		
Malt (n=8)	447.94	Average
	21.05	SD
	4.7	CV(%)
	99.54	Recovery (%)
Spike 1500 ppb		
Malt (n=8)	1556.1	Average
	46.32	SD
	2.97	CV(%)
	103.74	Recovery (%)
	101.64	Mean Recovery (%)

**Table 14. Recovery of Millet sample at two different levels.**

Spike 450 ppb		
Millet(n=8)	433.87	Average
	7.68	SD
	1.77	CV(%)
	96.4	Recovery (%)
Spike 1500 ppb		
Millet (n=8)	1458.55	Average
	61.23	SD
	4.19	CV(%)
	97.23	Recovery (%)
	96.81	Mean Recovery (%)

**Table 15. Recovery of Oats sample at two different levels.**

Spike 450 ppb		
Oats (n=8)	506.7	Average
	4.6	SD
	0.9	CV(%)
	112.6	Recovery (%)
Spike 1500 ppb		
Oats (n=8)	1582.1	Average
	88.92	SD
	5.6	CV(%)
	105.4	Recovery (%)
	110.5	Mean Recovery (%)

**Table 16. Recovery of Pasta sample at two different levels.**

Spike 450 ppb		
Pasta (n=8)	471	Average
	8	SD
	1.7	CV(%)
	104.66	Recovery (%)
Spike 1500 ppb		
Pasta (n=8)	1445.21	Average
	55.56	SD
	3.84	CV(%)
	96.34	Recovery (%)
	100.5	Mean Recovery (%)



**Table 17. Recovery of Pop corn sample at two different levels.**

Spike 450 ppb		
Pop Corn (n=8)	477	Average
	8.86	SD
	1.86	CV(%)
	106	Recovery (%)
Spike 1500 ppb		
Pop Corn (n=8)	1472	Average
	55	SD
	3.7	CV(%)
	98.13	Recovery (%)
	103.41	Mean Recovery (%)

**Table 18. Recovery of Rice Flour sample at two different levels.**

Spike 450 ppb		
Rice Flour (n=8)	406.9	Average
	14.5	SD
	3.56	CV(%)
	90.4	Recovery (%)
Spike 1500 ppb		
Rice Flour (n=8)	1437.2	Average
	45.82	SD
	3.18	CV(%)
	95.81	Recovery (%)
	93.1	Mean Recovery (%)

**Table 19. Recovery of Soybeans sample at two different levels.**

Spike 450 ppb		
Soybeans (n=8)	424.12	Average
	7.86	SD
	1.85	CV(%)
	94.25	Recovery (%)
Spike 1500 ppb		
Soybeans (n=8)	1445.14	Average
	55.35	SD
	38.13	CV(%)
	96.34	Recovery (%)
	95.3	Mean Recovery (%)

**Table 20. Recovery of Soybean Meal sample at two different levels.**

Spike 450 ppb		
Soybean Meal (n=8)	430.25	Average
	10.64	SD
	2.47	CV(%)
	95.6	Recovery (%)
Spike 1500 ppb		
Soybean Meal (n=8)	1486.2	Average
	87.22	SD
	5.86	CV(%)
	99.08	Recovery (%)
	97.34	Mean Recovery (%)

**Table 21. Recovery of Sunflower Meal sample at two different levels.**

Spike 450 ppb		
Sunflower Meal (n=8)	454.87	Average
	11.67	SD
	2.56	CV(%)
	101.1	Recovery (%)
Spike 1500 ppb		
Sunflower Meal (n=8)	1512.28	Average
	65.57	SD
	4.33	CV(%)
	100.8	Recovery (%)
	100.95	Mean Recovery (%)

**Table 22. Recovery of Wheat sample at two different levels.**

Spike 450 ppb		
Wheat (n=8)	434.6	Average
	9.64	SD
	2.22	CV(%)
	96.58	Recovery (%)
Spike 1500 ppb		
Wheat (n=8)	1488.3	Average
	54.2	SD
	3.64	CV(%)
	99.22	Recovery (%)
	97.9	Mean Recovery (%)

**Table 23. Recovery of Wheat Flour sample at two different levels.**

<b>Spike 450 ppb</b>		
<b>Wheat Flour (n=8)</b>	<b>475.75</b>	<b>Average</b>
	<b>9.86</b>	<b>SD</b>
	<b>2.07</b>	<b>CV(%)</b>
	<b>105.7</b>	<b>Recovery (%)</b>
<b>Spike 1500 ppb</b>		
<b>Wheat Flour (n=8)</b>	<b>1487.9</b>	<b>Average</b>
	<b>63.12</b>	<b>SD</b>
	<b>4.24</b>	<b>CV(%)</b>
	<b>99.19</b>	<b>Recovery (%)</b>
	<b>102.44</b>	<b>Mean Recovery (%)</b>

**Table 24. Recovery of White Rice sample at two different levels.**

<b>Spike 450 ppb</b>		
<b>White Rice (n=8)</b>	<b>468</b>	<b>Average</b>
	<b>6.44</b>	<b>SD</b>
	<b>1.37</b>	<b>CV(%)</b>
	<b>96.15</b>	<b>Recovery (%)</b>
<b>Spike 1500 ppb</b>		
<b>White Rice (n=8)</b>	<b>1462.1</b>	<b>Average</b>
	<b>57.8</b>	<b>SD</b>
	<b>3.95</b>	<b>CV(%)</b>
	<b>97.47</b>	<b>Recovery (%)</b>
	<b>96.81</b>	<b>Mean Recovery (%)</b>

**Table 25. Mean Recovery (%) of all ground matrices**

<b>Matrix</b>	<b>Mean Recovery (%)</b>
Barley	97.32
Beer Residue	101.26
Brow Rice	100.15
Buckwheat	100.03
Corn	100.56
Corn Flour	105.59
DDGS	99.93
Dried Brassica Integrifolia	90.63
Dried Gai Choy	91.47
Dried Palm	91.4
Malt	101.64
Millet	96.81
Oats	110.5
Pasta	100.5
Pop Corn	103.41
Rice Flour	93.1
Soybeans	95.3
Soybean Meal	97.34
Sunflower Meal	100.95
Wheat	97.9
Wheat Flour	102.44
White Rice	96.81
<b>MEAN</b>	<b>98.86</b>

### 3.3 Reproducibility

The coefficients of variation of reproducibility of the concentrations (ppb) (Table 26) of two different samples ran eight times in 8 different tests are reported:

**Table 26. Coefficients of Variation of the concentration (ppb) of two different samples ran in eight different tests.**

Sample (n=8)	Concentration (ppm)	
	MEAN	CV(%)
FAPAS MAIZE T04411QC	965.45	3.95
FAPAS MAIZE T04384QC	719.5	1.72

### 3.4 Performance Evaluation

#### i. Reference Materials


**Table 27. Recovery on samples prepared by FAPAS.**

Reference material	Lot number	Certified value (µg/kg)	Uncertainty (µg/kg)	Result (µg/kg)	Recovery (%)
FAPAS Maize T04411QC	1	962	310	965.45	100.36
FAPAS Maize T04384QC	87	752	251	719.5	104.5

## 4. References

- [1] Ntantasios AN, Arampatzis A, Voulgari D, Badra K, Papageorgiou G, Athanasiou SD and Gotsopoulos M. Innovative lateral flow method for the quantification of Aflatoxin M1. IDF DAIRY SUMMIT. 29 October-03 November 2017. Belfast. Northern Ireland. UK.
- [2] Papageorgiou G, Ntantasios AN, Voulgari D, Badra K, Gotsopoulos M and Athanasiou SD. An innovative symmetric lateral flow system for the quantification of Aflatoxin M1. 8th International Symposium on RAFA. 7-10 November 2017. Prague. Czech Republic.
- [3] M. Gkanas, Ch. Chatzoglou, K. Badra, Ch. Tsaridou, A.N. Ntantasios, G. Papageorgiou and S.D. Athanasiou. Uso di solventi non organici nell'analisi delle micotossine. Seminario AIA Laboratori e 20o ARAL SATA. 30-31 January 2018. Milan. Italy.
- [4] Drakouli S, Skliris A, Voulgari DL, Angeli E, Ntantasios AN, Papageorgiou G and Athanassiou SD, Estrazione unica in acqua, per la quantificazione di nove Micotossine usando la tecnologia Symmetric lateral flow. VI Congresso Nazionale: Micotossine e Tossine Vegetali nella filiera agro-alimentare, 10-12 June, 2019 Rome, Italy.
- [5] Tsaridou C, Badra K, Natsaridis N, Nikolopoulou E, Ntantasios AN, Papageorgiou G and Athanassiou SD, Estrazione unica in acqua, per la quantificazione di nove micotossine usando la tecnologia Bio-Shield Elisa. VI Congresso Nazionale: Micotossine e Tossine Vegetali nella filiera agro-alimentare, 10-12 June, 2019 Rome, Italy.
- [6] Drakouli S, Skliris A, Tziortziou M, Iliopoulou S, Natsaridis N, Papageorgiou G, Ntantasios AN and Athanassiou SD, Quantification of all Mycotoxins, using Symmetric lateral flow technology and one step multitoxin aqueous extraction. The World Mycotoxin Forum and the IUPAC International symposium on Mycotoxins, 14-16 October 2019, Belfast, Northern Ireland, UK.
- [7] Skliris A, Drakouli S, Tziortziou M, Voulgari DL, Iliopoulou S, Papageorgiou G, Zaralis K and Athanassiou SD, Symmetric lateral flow technology with one step Multitoxin aqueous extraction for the quantification of all Mycotoxins. 9th International Symposium on Recent Advances in Food Analysis, November 5-8, 2019, Prague, Czech Republic

Antonios Ntantasios  
Managing Director





[www.prognosis-biotech.com](http://www.prognosis-biotech.com)  
E: [info@prognosis-biotech.com](mailto:info@prognosis-biotech.com)

T: +30 2410 623922  
Farsalon 153 | 41335 Larissa, Greece